

ThermoSOLV-RX™

Protein Thermal Stability Kit

User Guide



Protocol: Identify Formulations that Maximize the Thermal Stability of a Protein Sample

Materials

Kit

- 1 tube TFluor™ Dye
- 1 Formulation Plate RX (96 x 170 μ L)

User-provided


- ca. 100 μ L purified protein solution (ca. 1 mg/mL)
- temperature scanning fluorescence plate reader such as a Q-PCR instrument (BioRad CFX, Roche LightCycler® 480, QuantStudio quantitative PCR)
- Q-PCR plate (preferably white)
- plate seal

Protocol

1. Dilute 2 μ L of TFluor Dye with 200 μ L of water and pipette 2 μ L of this solution onto the side of each well in a Q-PCR plate.
2. Pipette 1 μ L of protein sample into the opposite sides of all wells
3. Add 25 μ L of each formulation from the Formulation Plate RX. This combines the protein and the dye in each well.
4. Seal the Q-PCR plate and spin it (*i.e.* 5 min at 1000 rpm) to neatly collect all liquid in the center of the well.
5. Insert the plate into temperature scanning fluorescence plate reader. Run temperature scan (*i.e.* heating from 25°C to 90°C in 0.2°C steps equilibrating for 12 seconds for every step) while recording TFluor fluorescence at 570-580 nm (Excitation at 460-490nm).
6. Analyze the resulting fluorescence / temperature data and record the midpoint of thermal denaturation for each formulation. Compare thermal denaturation points and identify formulation that renders protein most temperature stable.


Variation: include known or putative small molecule ligands or co-factors (ATP, NAD/H, Zn²⁺ etc.) to protein buffer prior to analysis with the OptiTherm kit.

BioRad CFX is a trademark of Bio-Rad Laboratories, Inc.
QuantStudio is a trademark of Applied Biosystems
Roche LightCycler® is a trademark of Roche Diagnostics Corporation



Note & Troubleshooting

We advise to carry out a simple test prior to conducting the ThermoSOLV experiment to dial in the proper protein concentration and to identify a suitable detection range. This can be done by setting up a single well using the amounts suggested in this protocol (using any standard buffer). Consult the manual of the temperature scanning fluorescence plate reader to adjust the fluorescence emission signal to less than 20% of the maximal readout. Increase amounts of protein and dye if fluorescence signal is too low, decrease protein and dye amounts if fluorescence signal is too low.



ThermoSOLV-RX™

Product Information

Content:

- 1 x 96 well Formulation Plate RX
- 1 tube TFluor™ Dye
- Quick Start Guide and MSDS

Store at 4°C. Caution: TFluor™ Dye not yet fully tested (EU) WGK1 Combustible. Readily absorbed through skin. Target organ(s): Eyes, Skin. Hygroscopic. Product of USA. For R&D use only. Not for drug, household or other uses.

Purpose

ThermoSOLV-RX Protein Thermal Stability Kit

Systematic solution design and fluorescence-based stability assay for:

- **Thermal stabilization of protein samples**

For updated instructions and additional information please refer to www.stablebiologics.com

One ThermoSOLV-RX kit contains consumable materials to assay solution conditions for up to six (6) different protein samples.

Order Information

Order Cat #: SB-004-001

ThermoSOLV-RX Protein Thermal Stability Kit

Price: \$ 299 USD (3 pack discounted to \$750 USD)

Stable Biologics LLC
1500 1st Ave North
Birmingham, AL 35242
USA

Send order to: cweaver@stablebiologics.com

Reagent Listing

Well		Buffer#			Additive		Well		Buffer#			Additive					
#	Row Col	Conc	unit	pH	NAME	Conc	unit	#	Row Col	Conc	unit	pH	NAME	Conc	unit		
1	A 1				NaCl	60	mM	49	E 1				NaCl	60	mM		
2	A 2				Arg/Glu*	100	mM	50	E 2				Arg/Glu*	100	mM		
3	A 3				Arginine-HCl	30	mM	51	E 3				Arginine-HCl	30	mM		
4	A 4				Glycine	100	mM	52	E 4				Glycine	100	mM		
5	A 5				Poloxamer 188	0.2	%w/v	53	E 5				Poloxamer 188	0.2	%w/v		
6	A 6	Acetate	50	mM	5.0	EDTA	4	mM	54	E 6	Sodium phosphate	50	mM	6.5	EDTA	4	mM
7	A 7				Na bisulfate	6	mM	55	E 7				Na bisulfate	6	mM		
8	A 8				Sucrose	100	mM	56	E 8				Sucrose	100	mM		
9	A 9				Sorbitol	100	mM	57	E 9				Sorbitol	100	mM		
10	A 10				PEG 400	2	%w/v	58	E 10				PEG 400	2	%w/v		
11	A 11				Glycerol	6	%w/v	59	E 11				Glycerol	6	%w/v		
12	A 12	Ammonium sulfate	3	M				60	E 12	Sodium lactate	50	mM	6.5				
13	B 1				NaCl	60	mM	61	F 1				NaCl	60	mM		
14	B 2				Arg/Glu*	100	mM	62	F 2				Arg/Glu*	100	mM		
15	B 3				Arginine-HCl	30	mM	63	F 3				Arginine-HCl	30	mM		
16	B 4				Glycine	100	mM	64	F 4				Glycine	100	mM		
17	B 5				Poloxamer 188	0.2	%w/v	65	F 5				Poloxamer 188	0.2	%w/v		
18	B 6	Histidine	50	mM	6.0	EDTA	4	mM	66	F 6	Potassium phosphate	50	mM	7.0	EDTA	4	mM
19	B 7				Na bisulfate	6	mM	67	F 7				Na bisulfate	6	mM		
20	B 8				Sucrose	100	mM	68	F 8				Sucrose	100	mM		
21	B 9				Sorbitol	100	mM	69	F 9				Sorbitol	100	mM		
22	B 10				PEG 400	2	%w/v	70	F 10				PEG 400	2	%w/v		
23	B 11				Glycerol	6	%w/v	71	F 11				Glycerol	6	%w/v		
24	B 12	DMSO	5	%v/v				72	F 12	Na/K phosphate	50	mM	7.5	Tween 20	0.4	%w/v	
25	C 1				NaCl	60	mM	73	G 1				NaCl	60	mM		
26	C 2				Arg/Glu*	100	mM	74	G 2				Arg/Glu*	100	mM		
27	C 3				Arginine-HCl	30	mM	75	G 3				Arginine-HCl	30	mM		
28	C 4				Glycine	100	mM	76	G 4				Glycine	100	mM		
29	C 5				Poloxamer 188	0.2	%w/v	77	G 5				Poloxamer 188	0.2	%w/v		
30	C 6	Sodium Succinate	50	mM	6.0	EDTA	4	mM	78	G 6	Na/K phosphate	50	mM	7.5	EDTA	4	mM
31	C 7				Na bisulfate	6	mM	79	G 7				Na bisulfate	6	mM		
32	C 8				Sucrose	100	mM	80	G 8				Sucrose	100	mM		
33	C 9				Sorbitol	100	mM	81	G 9				Sorbitol	100	mM		
34	C 10				PEG 400	2	%w/v	82	G 10				PEG 400	2	%w/v		
35	C 11				Glycerol	6	%w/v	83	G 11				Glycerol	6	%w/v		
36	C 12	Original sample buffer						84	G 12				Benzyl alcohol	0.2	%w/v		
37	D 1				NaCl	60	mM	85	H 1				NaCl	60	mM		
38	D 2				Arg/Glu*	100	mM	86	H 2				Arg/Glu*	100	mM		
39	D 3				Arginine-HCl	30	mM	87	H 3				Arginine-HCl	30	mM		
40	D 4				Glycine	100	mM	88	H 4				Glycine	100	mM		
41	D 5				Poloxamer 188	0.2	%w/v	89	H 5				Poloxamer 188	0.2	%w/v		
42	D 6	Sodium citrate	50	mM	6.5	EDTA	4	mM	90	H 6	Tris	50	mM	7.5	EDTA	4	mM
43	D 7				Na bisulfate	6	mM	91	H 7				Na bisulfate	6	mM		
44	D 8				Sucrose	100	mM	92	H 8				Sucrose	100	mM		
45	D 9				Sorbitol	100	mM	93	H 9				Sorbitol	100	mM		
46	D 10				PEG 400	2	%w/v	94	H 10				PEG 400	2	%w/v		
47	D 11				Glycerol	6	%w/v	95	H 11				Glycerol	6	%w/v		
48	D 12	Glycine	50	mM	3.0	NaCl	500	mM	96	H 12	Tris	50	mM	8.5			

Note:

DMSO: dimethyl sulfoxide. EDTA: ethylenediaminetetraacetic acid. NaCl: sodium chloride. PEG: polyethylene glycol. # pH values for buffers used only; * Arginine-HCl/Glutamic acid, each amino acid is 50 mM
Well A12, B12 and C12 are for control experiments.